



medical **ECONET**
GERMANY



Only 5,4 kg

Wireless hand switch

Focal spot of 1,0 mm



mex+ AirRay

Ultra lightweight, portable X-ray Generator

- **HYBRID** – Device can be operated by internal battery or by external power supply
- Only 5,4kg for practical handling in mobile usage
- Equipped with high-performance battery (Lithium-Ion)
- Wireless charging or charging via USB connection
- Up to 60 exposures with only one battery charge
- Charging without removing the battery
- Exclusive remote control function by hand switch
- 10 programmable pre-set technique slots for kV and mAs
- Built-in ultra capacitor for steady power output

Specifications

Modell	meX+AirRay
Output	
Tube Vol. / Current	90kV / 20 mA
Voltage range	50 - 90 kV in 1 kV steps
mA Range	0,4mA - 20mA 25 steps
Max. output	1,6 kW at 80 kV
X-ray Tube	
Application	OX / 70 CEI
Focal Spot	1,0 mm
Collimator	
Max. Light field: 47 cm x 47 cm @ 100 cm SID	
Lamp	LED
Laser Pointer	dual laser
Features	
LCD-display, Hand switch remote control function	
Exposures full charge	max. 60 at 70 kV / 1 mAs
PROM Memory	10 APR
Weight	5,4 kg
Dimensions	325 x 268 x 250 mm

Standard Accessory

Wireless hand switch, hand switch and 2x rechargeable battery



Carrying Case




Display Features




Compact and easy to handle



Dimensions Generator (mm)



Dimensions Carrying Case (mm)



CE 0123

medical **ECONET**

Im Erlengrund 16, 46149 Oberhausen - Germany

t: +49 (0) 208 377 890-0 info@medical-econet.com
f: +49 (0) 208 377 890-55 www.medical-econet.com
facebook.com/medicaleconetgermany

Distributor

*Subject to change